DNA FRAGMENT TRANSLOCATION IN ARTIFICIAL SEA WATER THROUGH NANOPORES USING A PORTABLE MINI READER AND FLOW-CELL

David J. Niedzwiecki¹, Yung-Chien Chou², Michele Rossi³, Federico Thei³, Marija Drndic² ² Dep. Of Physics and Astronomy, University of Pennsylvania, Philadelphia, PA, USA ¹ Goeppert LLC, Philadelphia, PA, USA ³ Elements SRL, Cesena, Italy

LB109

Nanopore sensing is a powerful tool for the detection and characterization of biomolecules, such as nucleic acids. Solid-state nanopores act as single-molecule sensors that can function in a variety of harsh conditions, such as pH and temperature. The expected, resilient nature of solid-state nanopores makes them attractive potential candidates for taking this technology into the field to measure environmental samples for life detection applications and water quality monitoring. Towards this goal, here we measure artificial (synthetic) sea water samples spiked with DNA fragments at different concentrations to establish the limits of nanopore sensitivity in candidate environmental conditions. We use a compact, portable read-out device with miniaturized flow-cells, and obtain data with nanopores in 20-nm-thick, low-stress silicon nitride suspended on low-capacitance fused-silica (glass) chips at bandwidths up to 100 kHz. We measure 400 bp, 1000 bp, and 15000 bp DNA fragments, both mixed together and separately. We show that the fragments can be discriminated by analysis of current signal and dwell-time magnitudes.

Nanopore Mini Reader and Flow Cell



Abstract

- Low noise acquisition: ≈100 fA rms @ 1 KHz ≈380 fArms @ 10 KHz ≈3,2 pA rms @ 100
- 8 selectable sampling rates in the range 1.25 kHz 200 kHz;
- 4 input current rages: ± 200 pA (gain 2.25 G Ω), ± 2 nA (gain 225 M Ω); ± 20 nA (gain 22.5 M Ω);
- Voltage stimulus range between ±2000 mV, 1 mV minimum step;
- a fine automatic voltage offset compensation with 62,5 μ V step resolution;
- simple EDR software interface (also available for Mac OS) and data output format compatible with commercially available software (pClamp[®], Matlab[®], text);
- an online analysis module for real time analysis of data, as current histogram, Dwell time

Low-noise fused silica chips

- . Chip capacitance < 2 pF
- Fused silica substrate
- 20 nm LPCVD low-stress silicon nitride

. Low capacitance improves noise performance at high bandwidths [1]. Provided by Goeppert, LLC.

Corresponding author: David J. Niedzwiecki, Ph.D. - Goeppert LLC, Philadelphia, USA e-mail: dn@gppert.com

Measurements



Conclusions

- The combination of low-noise solid-state nanopores and the compact eNPR amplifier allowed for rapid characterization of DNA fragments. Analysis of DNA fragment events showed a separation of different fragment types in a Dwell-time vs. Current Amplitude scatter plot. Dwell-time of DNA fragments in the nanopore depended on DNA length.
- Spiking a sample of synthetic seawater with a DNA fragment showed that the seawater can be used as a detection medium without the need further buffering. Using different concentrations of a DNA spike demonstrates that detection in seawater is possible at picoMolar concentrations.

Acknowledgments

This work was funded in part by NASA SBIR Phase I Award 80NSSC18K1545 This work was carried out in part at the Singh Center for Nanotechnology, which is supported by the NSF National Nanotechnology Coordinated Infrastructure Program under grant NNCI-1542153.

Shekar, S.; Niedzwiecki, D. J.; Chien, C.-C.; Ong, P.; Fleischer, D. A.; Lin, J.; Rosenstein, J. K.; Drndić, M.; Shepard, K. L. Measurement of DNA Translocation Dynamics in a Solid-State Nanopore at 100 Ns Temporal Resolution. Nano Lett. 2016, 16 (7), 4483–4489, DOI: 10.1021/acs.nanolett.6b01661

Wanunu, M.; Morrison, W.; Rabin, Y.; Grosberg, A. Y.; Meller, A. Electrostatic focusing of unlabelled DNA into nanoscale pores using a salt gradient . Nat. Nanotechnol. 2010, 5, 160-165, DOI: 10.1038/nnano.2009.379



Visit us at booth 400

References