
EDA Installation and Licensing

Revision History

Date	Version	Description
05/11/2019	1.0	First version of document



Installation

EDA is part of the EDR3 installation bundle, so in order to get EDA just download EDR3 software from our [website](#) and launch the installation program.


Licensing

What's included in the license

Feature	Not licensed	Licensed
Import/export dat, abf and csv	V	V
Data and tags navigation	V	V
Digital low pass filters	X	V
Offline analysis	X	V

Activating the license

Launch EDA, click the question mark and select the License option





In the windows that will open copy the serial code (you can use the Copy button to be sure to copy all of the characters in the code) and send it to info@elements-ic.com together with your device ID (S/N).

Elements will send you the Activation Code to be copied in the lower field.

Activate EDA

To receive the activation code for EDA please send the serial code shown below to info@elements-ic.com (in the email add also your device ID).
The Elements team will reply to you as soon as possible providing the activation code.

Serial Code:

Activation Code:

Please, notice that the activation code is machine locked, so if you want to use EDA on multiple computers you will have to ask for the activation code for each machine.