




e16HC (built-in cage) dimensions

PCB details



PCB pinout

pin1 = reference	pin6 = in 12	pin11 = in 7	pin16 = in 2
pin2 = in 16	pin7 = in 11	pin12 = in 6	pin17 = in 1
pin3 = in 15	pin8 = in 10	pin13 = in 5	pin18 = reference
pin4 = in 14	pin9 = in 9	pin14 = in 4	pin23 = ground
pin5 = in 13	pin10 = in 8	pin15 = in 3	

Maximum available volume

